

WHAT IS THE EYEBAG?

When trying to unblock the Meibomian glands, what matters is that the heating is gentle and continuous and the EyeBag does just that. Conventional warm compresses are either too hot to start with and cool down very quickly. If you think about when you melt chocolate, it needs to be done gently and it can't get too hot otherwise the consistency gets too lumpy!

The MGDRx EyeBag is a silk and suedex device which can be re-heated up to 200 times in a microwave. Simply heat for 30 or 40 seconds on full power and the EyeBag will stay warm for up to 10 minutes.

The eye bag was designed by an Ophthalmologist and it works by gently heating up the matter in the blocked glands to the right temperature making unblocking the glands much more easy.

Conditions which usually benefit from regular twice daily EyeBag treatment are:

- Meibomian Gland Dysfunction (MGD)
- Blepharitis
- Dry Eye Syndrome
- Styes
- Eyelid Cysts
- Contact Lens Dry Eye
- Corneal Laser Dry Eye
- Rosacea



Symptoms which are relieved by EyeBag treatment are:

- Grittiness
- Tiredness
- Redness
- Blurry Vision
- Burning Discomfort
- Foreign body sensation
- Dry eye
- Eye irritation

